

De-Blister Machine (DB500)

PRODUCT CODE: 100-62M

Our extensive range of products has been created by pharmacists for pharmacists and we pride ourselves on our individual service and personal touch.

You might also be interested in:

- Our range of SureMed by Omnicell multi-dose and single-dose pill packs.
- The MTS 350 filling and sealing machine

“The quickest and easiest way to de-blister medication.”

Take the strain and pain out of removing tablets from their packaging with our de-blister machine (DB500), designed specifically to transform the time-consuming task of de-blistering.

This compact desktop machine not only reduces the risk of repetitive injuries and sore fingers from de-blistering by hand, it's also quick to set-up for different medication packs with no tools needed.

Designed to work with the vast majority of medication packs used in today's pharmacy, the DB500 can de-blister up to 280 capsules or tablets per minute.

It's simple to clean and maintain and comes with a safety clamp to secure the unit to your workstation when in use.

We'll provide a comprehensive DVD and instruction manual to help you use the DB500 with confidence. It's so straightforward to use that operators can be trained in a matter of minutes.

Like all Omnicell products, the DB500 is designed to support you in providing a valuable service to your customers and to support them in taking the correct medication.

Want to find out more?

We are the world's leading provider of medication adherence packaging, solutions and technology so contact us today using the number below to find out what we can do for you.



The DB500 transforms the time consuming task of de-blistering.



The DB500 can de-blister up to 280 capsules or tablets per minute

Omnicell[®]

Contact us today:

Call: +44 (0)161 413 5333

Fax: +44 (0)161 775 6658

Email: ordersuk@omnicell.com

www.omnicell.co.uk